

LGALS2

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:50 - 1:200

Calculated MW: 15kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LGALS2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HL14;

Catalog #: A2734

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3957

Isotype: IgG

Swiss Prot: P05162

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a soluble beta-galactoside binding lectin. The encoded protein is found as a homodimer and can bind to lymphotoxin-alpha. A single nucleotide polymorphism in an intron of this gene can alter the transcriptional level of the protein, with a resultant increased risk of myocardial infarction.

To place an order, please [Click HERE](#).