

GAS7

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:100

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human GAS7

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MLL/GAS7;

Catalog #: A2736

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8522

Isotype: IgG

Swiss Prot: O60861

Purity: Affinity purification

For research use only.

Background:

Growth arrest-specific 7 is expressed primarily in terminally differentiated brain cells and predominantly in mature cerebellar Purkinje neurons. GAS7 plays a putative role in neuronal development. Several transcript variants encoding proteins which vary in the N-terminus have been described.

To place an order, please [Click HERE](#).