

GDF3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:25 - 1:100

Calculated MW: 41kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GDF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KFS3; MCOP7; MCOPCB6;

Catalog #: A2738

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9573

Isotype: IgG

Swiss Prot: Q9NR23

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the bone morphogenetic protein (BMP) family and the TGF-beta superfamily. This group of proteins is characterized by a polybasic proteolytic processing site which is cleaved to produce a mature protein containing seven conserved cysteine residues. The members of this family are regulators of cell growth and differentiation in both embryonic and adult tissues.

To place an order, please [Click HERE](#).