

GLP2R

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:63kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GLP2R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A2740

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9340

Isotype:IgG

Swiss Prot:O95838

Purity:Affinity purification

For research use only.

Background:

This gene encodes a G protein-coupled receptor that is closely related to the glucagon receptor and binds to glucagon-like peptide-2 (GLP2). Signalling through GLP2 stimulates intestinal growth and increases villus height in the small intestine, concomitant with increased crypt cell proliferation and decreased enterocyte apoptosis.

To place an order, please [Click HERE](#).