

GPC6

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 63kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GPC6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OMIMD1;

Catalog #: A2741

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10082

Isotype: IgG

Swiss Prot: Q9Y625

Purity: Affinity purification

For research use only.

Background:

The glypicans comprise a family of glycosylphosphatidylinositol-anchored heparan sulfate proteoglycans, and they have been implicated in the control of cell growth and cell division. The glypican encoded by this gene is a putative cell surface coreceptor for growth factors, extracellular matrix proteins, proteases and anti-proteases.

To place an order, please [Click HERE](#).