

GRB14

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:100 IHC 1:25 - 1:100

Calculated MW:61kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GRB14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A2742

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2888

Isotype:IgG

Swiss Prot:Q14449

Purity:Affinity purification

For research use only.

Background:

The product of this gene belongs to a small family of adapter proteins that are known to interact with a number of receptor tyrosine kinases and signaling molecules. This gene encodes a growth factor receptor-binding protein that interacts with insulin receptors and insulin-like growth-factor receptors. This protein likely has an inhibitory effect on receptor tyrosine kinase signaling and, in particular, on insulin receptor signaling. This gene may play a role in signaling pathways that regulate growth and metabolism. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).