

HMGB4

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 22kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HMGB4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

dJ1007G16.5;

To place an order, please [Click HERE](#).

Catalog #: A2744

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 127540

Isotype: IgG

Swiss Prot: Q8WW32

Purity: Affinity purification

For research use only.