

TRAP1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:80kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human TRAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

HSP75; HSP 75; HSP90L; TRAP-1;

Catalog #:A2748

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10131

Isotype:IgG

Swiss Prot:Q12931

Purity:Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial chaperone protein that is member of the heat shock protein 90 (HSP90) family. The encoded protein has ATPase activity and interacts with tumor necrosis factor type I. This protein may function in regulating cellular stress responses. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).