

KLK15

Reactivity: Human Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KLK15

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ACO; HSRNASPH;

Catalog #: A2757

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55554

Isotype: IgG

Swiss Prot: Q9H2R5

Purity: Affinity purification

For research use only.

Background:

Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. In prostate cancer, this gene has increased expression, which indicates its possible use as a diagnostic or prognostic marker for prostate cancer. The gene contains multiple polyadenylation sites and alternative splicing results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).