

MAGEB4

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:15 - 1:50

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MAGEB4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CT3.6;

Catalog #: A2760

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4115

Isotype: IgG

Swiss Prot: O15481

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the MAGEB gene family. The members of this family have their entire coding sequences located in the last exon, and the encoded proteins show 50 to 68% sequence identity to each other. The promoters and first exons of the MAGEB genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEB genes are clustered on chromosome Xp22-p21. This gene sequence ends in the first intron of MAGEB1, another family member. This gene is expressed in testis.

To place an order, please [Click HERE](#).