

MATN4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 68kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MATN4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A2761

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8785

Isotype: IgG

Swiss Prot: O95460

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of von Willebrand factor A domain-containing protein family. The proteins of this family are thought to be involved in the formation of filamentous networks in the extracellular matrices of various tissues. This family member is thought to play a role in reorganizing and regenerating the corneal matrix in granular and lattice type I dystrophies. It may also be involved in wound healing in the dentin-pulp complex. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).