

MMP28

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MMP28

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MM28; MMP25; MMP-25; MMP-28; EPILYSIN;

Catalog #: A2764

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79148

Isotype: IgG

Swiss Prot: Q9H239

Purity: Affinity purification

For research use only.

Background:

Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix for both normal physiological processes, such as embryonic development, reproduction and tissue remodeling, and disease processes, such as asthma and metastasis. This gene encodes a secreted enzyme that degrades casein. Its expression pattern suggests that it plays a role in tissue homeostasis and in wound repair. Alternative splicing of this gene results in multiple transcript variants.

To place an order, please [Click HERE](#).