

MYOZ1

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MYOZ1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CS-2; FATZ; MYOZ;

Catalog #: A2767

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 58529

Isotype: IgG

Swiss Prot: Q9NP98

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is primarily expressed in the skeletal muscle, and belongs to the myozenin family. Members of this family function as calcineurin-interacting proteins that help tether calcineurin to the sarcomere of cardiac and skeletal muscle. They play an important role in modulation of calcineurin signaling.

To place an order, please [Click HERE](#).