

NAP1L1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:45kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NAP1L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

NRP; NAP1; NAP1L;

Catalog #:A2769

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4673

Isotype:IgG

Swiss Prot:P55209

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the nucleosome assembly protein (NAP) family. This protein participates in DNA replication and may play a role in modulating chromatin formation and contribute to the regulation of cell proliferation. Alternative splicing of this gene results in several transcript variants; however, not all have been fully described.

To place an order, please [Click HERE](#).