

MAP3K14

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:25 - 1:100

Calculated MW: 104kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MAP3K14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HS; NIK; HSNIK; FTDCR1B;

Catalog #: A2773

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9020

Isotype: IgG

Swiss Prot: Q99558

Purity: Affinity purification

For research use only.

Background:

This gene encodes mitogen-activated protein kinase kinase kinase 14, which is a serine/threonine protein-kinase. This kinase binds to TRAF2 and stimulates NF-kappaB activity. It shares sequence similarity with several other MAPKK kinases. It participates in an NF-kappaB-inducing signalling cascade common to receptors of the tumour-necrosis/nerve-growth factor (TNF/NGF) family and to the interleukin-1 type-I receptor.

To place an order, please [Click HERE](#).