

POU2F2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human POU2F2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

OCT2; OTF2; Oct-2;

Catalog #: A2776

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5452

Isotype: IgG

Swiss Prot: P09086

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a homeobox-containing transcription factor of the POU domain family. The encoded protein binds the octamer sequence 5'-ATTGTCAT-3', a common transcription factor binding site in immunoglobulin gene promoters. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).