

OPCML

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 38kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human OPCML

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OPCM; OBCAM; IGLON1;

Catalog #: A2778

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4978

Isotype: IgG

Swiss Prot: Q14982

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the IgLON subfamily in the immunoglobulin protein superfamily. The encoded protein is localized in the plasma membrane and may have an accessory role in opioid receptor function. This gene has an ortholog in rat and bovine. The opioid binding-cell adhesion molecule encoded by the rat gene binds opioid alkaloids in the presence of acidic lipids, exhibits selectivity for mu ligands and acts as a GPI-anchored protein. Since the encoded protein is highly conserved in species during evolution, it may have a fundamental role in mammalian systems. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).