

## FAM3A

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:1000 IHC 1:50 - 1:200

**Calculated MW:**25kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human FAM3A

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

DLD; 2.19; XAP-7; DXS560S;

**Catalog #:**A2784

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**60343

**Isotype:**IgG

**Swiss Prot:**P98173

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a cytokine-like protein. The expression of this gene may be regulated by peroxisome proliferator-activated receptor gamma, and the encoded protein may be involved in the regulation of glucose and lipid metabolism. Alternative splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*