

APAF1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human APAF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bi

Synonym:

CED4; APAF-1

Catalog #: A2800

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 317

Isotype: IgG

Swiss Prot: O14727

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cytoplasmic protein that initiates apoptosis. This protein contains several copies of the WD-40 domain, a caspase recruitment domain (CARD), and an ATPase domain (NB-ARC). Upon binding cytochrome c and dATP, this protein forms an oligomeric apoptosome. The apoptosome binds and cleaves caspase 9 preproprotein, releasing its mature, activated form. Activated caspase 9 stimulates the subsequent caspase cascade that commits the cell to apoptosis. Alternative splicing results in several transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).