

RAB38

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human RAB38

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAB38, Ras related protein Rab 38, Melanoma antigen NY MEL 1, NY MEL 1, rrGTPbp;

Catalog #: A2813

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23682

Isotype: IgG

Swiss Prot: P57729

Purity: Affinity purification

For research use only.

Background:

RAB38 (RAB38, member RAS oncogene family) is a protein-coding gene. Diseases associated with RAB38 include melanoma. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is RABL2B. May be involved in melanosomal transport and docking. Involved in the proper sorting of TYRP1. Plays a role in the maturation of phagosomes that engulf pathogens, such as S.aureus and M.tuberculosis.

To place an order, please [Click HERE](#).