

RAB3C

Reactivity: Human Mouse

Tested applications: WB IHC IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:10 - 1:50 IP 1:20 - 1:50

Calculated MW: 26kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human RAB3C

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A2823

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 115827

Isotype: IgG

Swiss Prot: Q96E17

Purity: Affinity purification

For research use only.

Background:

RAB proteins, such as RAB3C, are small GTPases implicated in the regulation of vesicular trafficking between membrane compartments (Cheng et al., 2002)

To place an order, please [Click HERE](#).