

ARL2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human ARL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

ARFL2;

Catalog #: A2830

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 402

Isotype: IgG

Swiss Prot: P36404

Purity: Affinity purification

For research use only.

Background:

This gene encodes a small GTP-binding protein of the RAS superfamily which functions as an ADP-ribosylation factor (ARF). The encoded protein is one of a functionally distinct group of ARF-like genes.

To place an order, please [Click HERE](#).