

HAVCR1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human HAVCR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TIM; KIM1; TIM1; CD365; HAVCR; KIM-1; TIM-1; TIMD1; TIMD-1; HAVCR-1;

Catalog #: A2831

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26762

Isotype: IgG

Swiss Prot: Q96D42

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a membrane receptor for both human hepatitis A virus (HHAV) and TIMD4. The encoded protein may be involved in the moderation of asthma and allergic diseases. The reference genome represents an allele that retains a MTTVP amino acid segment that confers protection against atopy in HHAV seropositive individuals. Alternative splicing of this gene results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 4, 12 and 19.

To place an order, please [Click HERE](#).