

BDKRB2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 44kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BDKRB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

B2R; BK2; BK-2; BKR2; BRB2;

Catalog #: A2844

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 624

Isotype: IgG

Swiss Prot: P30411

Purity: Affinity purification

For research use only.

Background:

This gene encodes a receptor for bradykinin. The 9 aa bradykinin peptide elicits many responses including vasodilation, edema, smooth muscle spasm and pain fiber stimulation. This receptor associates with G proteins that stimulate a phosphatidylinositol-calcium second messenger system. Alternate start codons result in two isoforms of the protein.

To place an order, please [Click HERE](#).