

BCL2

Reactivity: Human Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human BCL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bi

Synonym:

Bcl-2;PPP1R50;

Catalog #: A2845

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 596

Isotype: IgG

Swiss Prot: P10415

Purity: Affinity purification

For research use only.

Background:

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Two transcript variants, produced by alternate splicing, differ in their C-terminal ends.

To place an order, please [Click HERE](#).