

CDC37

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:44kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CDC37

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

P50CDC37;

Catalog #:A2859

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11140

Isotype:IgG

Swiss Prot:Q16543

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is highly similar to Cdc 37, a cell division cycle control protein of *Saccharomyces cerevisiae*. This protein is a molecular chaperone with specific function in cell signal transduction. It has been shown to form complex with Hsp90 and a variety of protein kinases including CDK4, CDK6, SRC, RAF-1, MOK, as well as eIF2 alpha kinases. It is thought to play a critical role in directing Hsp90 to its target kinases.

To place an order, please [Click HERE](#).