

CEBPD

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CEBPD

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CELF; CRP3; C/EBP-delta; NF-IL6-beta;

Catalog #: A2862

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1052

Isotype: IgG

Swiss Prot: P49716

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this intronless gene is a bZIP transcription factor which can bind as a homodimer to certain DNA regulatory regions. It can also form heterodimers with the related protein CEBP-alpha. The encoded protein is important in the regulation of genes involved in immune and inflammatory responses, and may be involved in the regulation of genes associated with activation and/or differentiation of macrophages. The cytogenetic location of this locus has been reported as both 8p11 and 8q11.

To place an order, please [Click HERE](#).