

CHGB

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:100 - 1:500

Calculated MW:78kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CHGB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SCG1;

Catalog #:A2864

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1114

Isotype:IgG

Swiss Prot:P05060

Purity:Affinity purification

For research use only.

Background:

This gene encodes a tyrosine-sulfated secretory protein abundant in peptidergic endocrine cells and neurons. This protein may serve as a precursor for regulatory peptides.

To place an order, please [Click HERE](#).