

CSNK1A1

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:20 - 1:100

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CSNK1A1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CK1; CK1a; CK1α; HLCDGP1; PRO2975; HEL-S-77p;

To place an order, please [Click HERE](#).

Catalog #: A2868

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1452

Isotype: IgG

Swiss Prot: P48729

Purity: Affinity purification

For research use only.