

DPYSL3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 62kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human DPYSL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DRP3; ULIP; CRMP4; DRP-3; LCRMP; CRMP-4; ULIP-1;

Catalog #: A2877

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1809

Isotype: IgG

Swiss Prot: Q14195

Purity: Affinity purification

For research use only.

Background:

Collapsin response mediator protein 4 is an enzyme that in humans is encoded by the CRMP4 gene. Mostly expressed in heart; skeletal muscle; fetal brain and spinal cord. It has phosphorylation sites on Ser-522 by DYRK2 promotes subsequent phosphorylation on Thr-509; Thr-514 and Ser-518 by GSK3.

To place an order, please [Click HERE](#).