

DPYSL5

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 61kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human DPYSL5

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRAM; CRMP5; Ulip6; CRMP-5;

Catalog #: A2878

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 56896

Isotype: IgG

Swiss Prot: Q9BPU6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the CRMP (collapsing response mediator protein) family thought to be involved in neural development. Antibodies to the encoded protein were found in some patients with neurologic symptoms who had paraneoplastic syndrome. A pseudogene of this gene is found on chromosome 11. Multiple alternatively spliced variants, encoding the same protein, have been identified.

To place an order, please [Click HERE](#).