

SKI

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1000

Calculated MW: 80kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SKI

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

c-SKI;

Catalog #: A2879

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6497

Isotype: IgG

Swiss Prot: P12755

Purity: Affinity purification

For research use only.

Background:

The Ski protein is a nuclear protooncoprotein that is associated with tumors at high cellular concentrations. It has been shown to interfere with normal cellular functioning by both directly impeding expression of certain genes inside the nucleus of the cell and disrupting signaling proteins that activate genes. This protein may play a role in terminal differentiation of skeletal muscle cells but not in the determination of cells to the myogenic lineage. May functions in TGF-beta signaling.

To place an order, please [Click HERE](#).