

CUL2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 87kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CUL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

To place an order, please [Click HERE](#).

Catalog #: A2881

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8453

Isotype: IgG

Swiss Prot: Q13617

Purity: Affinity purification

For research use only.