

GJA9

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human GJA9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CX58; CX59; GJA10;

Catalog #: A2884

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81025

Isotype: IgG

Swiss Prot: P57773

Purity: Affinity purification

For research use only.

Background:

Connexins, such as GJA9, are involved in the formation of gap junctions, intercellular conduits that directly connect the cytoplasm of contacting cells. Each gap junction channel is formed by docking of 2 hemichannels, each of which contains 6 connexin subunits (Sohl et al., 2003)

To place an order, please [Click HERE](#).