

AQP10

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human AQP10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AQPA_HUMAN;

Catalog #: A2888

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 89872

Isotype: IgG

Swiss Prot: Q96PS8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the aquaglyceroporin family of integral membrane proteins. Members of this family function as water-permeable channels in the epithelia of organs that absorb and excrete water. This protein was shown to function as a water-selective channel, and could also permeate neutral solutes such as glycerol and urea.

To place an order, please [Click HERE](#).