

CX3CR1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 44kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CX3CR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

V28l;CCRL1;GPR13;CMKDR1;GPRV28;CMKBRL1;

Catalog #: A2890

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1524

Isotype: IgG

Swiss Prot: P49238

Purity: Affinity purification

For research use only.

Background:

Fractalkine is a transmembrane protein and chemokine involved in the adhesion and migration of leukocytes. The protein encoded by this gene is a receptor for fractalkine. The encoded protein also is a coreceptor for HIV-1, and some variations in this gene lead to increased susceptibility to HIV-1 infection and rapid progression to AIDS. Four transcript variants encoding two different isoforms have been found for this gene.

To place an order, please [Click HERE](#).