

EDNRA

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human EDNRA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ETA; ET-A; ETAR; ETRA; MFDA; ETA-R; hET-AR;

Catalog #: A2907

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1909

Isotype: IgG

Swiss Prot: P25101

Purity: Affinity purification

For research use only.

Background:

This gene encodes the receptor for endothelin-1, a peptide that plays a role in potent and long-lasting vasoconstriction. This receptor associates with guanine-nucleotide-binding (G) proteins, and this coupling activates a phosphatidylinositol-calcium second messenger system. Polymorphisms in this gene have been linked to migraine headache resistance. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).