

EDNRB

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 50kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human EDNRB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ETB;ET-B;ETBR;ETRB;HSCR;WS4A;ABCDS;ET-BR;HSCR2;

Catalog #: A2908

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1910

Isotype: IgG

Swiss Prot: P24530

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a G protein-coupled receptor which activates a phosphatidylinositol-calcium second messenger system. Its ligand, endothelin, consists of a family of three potent vasoactive peptides: ET1, ET2, and ET3. Studies suggest that the multigenic disorder, Hirschsprung disease type 2, is due to mutations in the endothelin receptor type B gene. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).