

SFRP1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:100 - 1:500

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SFRP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FRP; FRP1; FrzA; FRP-1; SARP2;

Catalog #: A2911

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6422

Isotype: IgG

Swiss Prot: Q8N474

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. Members of this family act as soluble modulators of Wnt signaling; epigenetic silencing of SFRP genes leads to deregulated activation of the Wnt-pathway which is associated with cancer. This gene may also be involved in determining the polarity of photoreceptor cells in the retina.

To place an order, please [Click HERE](#).