

FOXD3

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 48kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human FOXD3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AIS1; HFH2; VAMAS2; Genesis;

Catalog #: A2926

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27022

Isotype: IgG

Swiss Prot: Q9UJU5

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.

To place an order, please [Click HERE](#).