

FOXL1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 36kD

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human FOXL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

FKH6;FKHL11;FREAC7;

Catalog #: A2931

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2300

Isotype: IgG

Swiss Prot: Q12952

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the forkhead/winged helix-box (FOX) family of transcription factors. FOX transcription factors are characterized by a distinct DNA-binding forkhead domain and play critical roles in the regulation of multiple processes including metabolism, cell proliferation and gene expression during ontogenesis.

To place an order, please [Click HERE](#).