

FOXL2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human FOXL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BPES; PFRK; POF3; BPES1; PINTO;

Catalog #: A2932

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 668

Isotype: IgG

Swiss Prot: P58012

Purity: Affinity purification

For research use only.

Background:

This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-binding domain and may play a role in ovarian development and function. Mutations in this gene are a cause of blepharophimosis syndrome and premature ovarian failure 3.

To place an order, please [Click HERE](#).