

FOXP1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 75kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human FOXP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MFH; QRF1; 12CC4; hFKH1B; HSPC215;

Catalog #: A2935

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27086

Isotype: IgG

Swiss Prot: Q9H334

Purity: Affinity purification

For research use only.

Background:

This gene belongs to subfamily P of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in the regulation of tissue- and cell type-specific gene transcription during both development and adulthood. Forkhead box P1 protein contains both DNA-binding- and protein-protein binding-domains. This gene may act as a tumor suppressor as it is lost in several tumor types and maps to a chromosomal region (3p14.1) reported to contain a tumor suppressor gene(s). Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).