

CCKBR

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 48kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CCKBR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GASR; CCK-B; CCK2R;

Catalog #: A2941

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 887

Isotype: IgG

Swiss Prot: P51959

Purity: Affinity purification

For research use only.

Background:

This gene encodes a G-protein coupled receptor for gastrin and cholecystokinin (CCK), regulatory peptides of the brain and gastrointestinal tract. This protein is a type B gastrin receptor, which has a high affinity for both sulfated and nonsulfated CCK analogs and is found principally in the central nervous system and the gastrointestinal tract. A misspliced transcript variant including an intron has been observed in cells from colorectal and pancreatic tumors.

To place an order, please [Click HERE](#).