

CMTM5

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CMTM5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CKLFSF5;

Catalog #: A2942

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 116173

Isotype: IgG

Swiss Prot: Q96DZ9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the chemokine-like factor superfamily. This family of genes encodes multi-pass membrane proteins that are similar to both the chemokine and the transmembrane 4 superfamilies of signaling molecules. The encoded protein may exhibit tumor suppressor activity. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).