

CMTM3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 20kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CMTM3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BNAS2; CKLFSF3

Catalog #: A2943

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 123920

Isotype: IgG

Swiss Prot: Q96MX0

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies of signaling molecules. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 16. Alternatively spliced transcript variants containing different 5' UTRs, but encoding the same protein, have been identified.

To place an order, please [Click HERE](#).