

GSN

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 90kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GSN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADF; AGEL

Catalog #: A2952

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2934

Isotype: IgG

Swiss Prot: P06396

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene binds to the "plus" ends of actin monomers and filaments to prevent monomer exchange. The encoded calcium-regulated protein functions in both assembly and disassembly of actin filaments. Defects in this gene are a cause of familial amyloidosis Finnish type (FAF). Multiple transcript variants encoding several different isoforms have been found for this gene.

To place an order, please [Click HERE](#).