

GFRA2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GFRA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NTNRA; RETL2; TRNR2; GDNFRB; NRTNR-ALPHA;

Catalog #: A2954

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2675

Isotype: IgG

Swiss Prot: O00451

Purity: Affinity purification

For research use only.

Background:

Glial cell line-derived neurotrophic factor (GDNF) and neurturin (NTN) are two structurally related, potent neurotrophic factors that play key roles in the control of neuron survival and differentiation. The protein encoded by this gene is a member of the GDNF receptor family. It is a glycosylphosphatidylinositol(GPI)-linked cell surface receptor for both GDNF and NTN, and mediates activation of the RET tyrosine kinase receptor. This encoded protein acts preferentially as a receptor for NTN compared to its other family member, GDNF family receptor alpha 1. This gene is a candidate gene for RET-associated diseases. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).