

GFRA3

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:45kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GFRA3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GDNFR3;

Catalog #:A2955

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2676

Isotype:IgG

Swiss Prot:O60609

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a glycosylphosphatidylinositol(GPI)-linked cell surface receptor and a member of the GDNF receptor family. It forms a signaling receptor complex with RET tyrosine kinase receptor and binds the ligand, artemin.

To place an order, please [Click HERE](#).