

PROKR2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 44kD

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human PROKR2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HH3; KAL3; PKR2; GPRg2; GPR73b; GPR73L1; dJ680N4.3;

Catalog #: A2960

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 128674

Isotype: IgG

Swiss Prot: Q8NFJ6

Purity: Affinity purification

For research use only.

Background:

Prokineticins are secreted proteins that can promote angiogenesis and induce strong gastrointestinal smooth muscle contraction. The protein encoded by this gene is an integral membrane protein and G protein-coupled receptor for prokineticins. The encoded protein is similar in sequence to GPR73, another G protein-coupled receptor for prokineticins.

To place an order, please [Click HERE](#).